



专注牙科设备常用配件研发制造
Focus on R & D and manufacturing
of common dental equipment accessories



产品型录 PRODUCT CATALOG


Suzhou ABC Electronics Technology Co.,Ltd





Focus on R & D and manufacturing of common dental equipment accessories

我们的产品
OUR PRODUCT

	01	负压控制阀 Vacuum Suction Valve		09	管夹 Pipe Clamp
	02	选位阀 Place Selection Valve		09	管箍 Shroud Ring
	06	控制电位器 Potentiometer		10	管套 Pipe Shroud
	06	牙科抽吸管 Dental Suction Pipes		11	接嘴 Nipple
	07	直接头 Straight Connector		11	卡扣 U-Clip
	08	内接头 Internal Connector		11	O型密封圈 O Seal Ring



12

接头套装
Connectors Kit



14

接头套装
Connectors Kit



14

强力吸唾头
High Volume Ejector



14

可调式吸唾头
Adjustable Suction
Connector



15

吸唾接头
Suction Connector



16

可调式吸唾头
Adjustable Suction
Connector



17

可调式吸唾头
Adjustable Suction
Connector



18

吸唾转换头
Suction Adapter



19

吸管接头
Suction Connector



20

直通接头
Straight Joint



20

45°/90°接头
Connector

我们的产品
OUR PRODUCT



20

转接头
Adapter



20

排气接头
Air Outlet Connector




21

单向阀
Check Valve



21

过滤网接头
Filter Net Connector



22

三通接头
T-Junction



22

封盖
Seal Cap



23

多通接头
Multi-Connectors



23

抽吸释压阀
Vacuum Release Valve



23

截流开关
Air Flow Adjust Valve



24

装饰管套
Decoration Casing Pipe



24

回气瓶
Air Return Bottle



24

排污管接头
Drain Pipe Connector

③

专注牙科设备常用配件研发制造

Focus on R & D and manufacturing of common dental equipment accessories



25

水过滤器
Water Filter




25

水过滤网
Water Filter Net



26

过滤筒
Filter Cartridge



26

过滤器面板
Filter Panel



27

过滤器接头
Filter Connector




29

过滤网
Filter Net



29

过滤网杆
Disposable Filter Rod



30

过滤网组合套
Disposable Filter Net



31

过滤器组合套
Filter Sets



33

过滤器组合套
Filter Sets



35

过滤器组合套
Filter Sets



36

过滤器组合套
Filter Sets

产品型录
PRODUCT CATALOG


单位UNIT: mm

负压控制阀 Vacuum Suction Valve



J38-0139

带控制电位器负压控制阀
Vacuum Suction Valve with Potentiometer



J38-001

电控式负压控制阀
Electronic Vacuum Suction Valve



J38-010

气控式负压控制阀
Pneumatic Vacuum Suction Valve

	J38-0139	J38-001	J38-010
延时功能 Time-delay Protection	●	○	○
无负压启动 Start without vacuum pressure	●	●	●
全面隔绝 Complete Isolation	●	●	●
安装灵便 Easy Assembly	●	●	●
快速装卸 Quick Dis-assembly	●	●	●
电源输入 Power Input	●	●	○
抗扰设计 Anti-interference	●	○	○

*说明：●是 ○否，以上产品配套齐全，适用于干抽式或半湿抽式的中央负压抽吸设备。
Note：● Yes ○ No, Complete Sets for Assembly, VALVE shall be suitable for dry suction unit and semi-dry suction unie.

产品型录
PRODUCT CATALOG

单位UNIT: mm

选位阀Place Selection Valve




J06-0107

带控制电位器选位阀
Place Selection Valve With Potentiometer



J06-001

电控式选位阀
Electronic Place Selection Valve



J25-001

气控式选位阀
Pneumatic Type Place Selection Valve

	J06-0107	J06-001	J25-001
延时功能 Time-delay Protection	●	○	○
无负压启动 Start without vacuum pressure	○	○	○
全面隔绝 Complete Isolation	○	○	○
安装灵便 Easy Assembly	○	○	○
快速装卸 Quick Dis-assembly	○	○	○
电源输入 Power Input	●	●	○
抗扰设计 Anti-interference	●	○	○

*说明: ●是 ○否, 以上产品配套齐全, 适用于干抽式或半湿抽式的中央负压抽吸设备。
Note: ● Yes ○ No, Complete Sets for Assembly, VALVE shall be suitable for dry suction unit and semi-dry suction unie.

产品型录
PRODUCT CATALOG

单位UNIT: mm

产品优势说明Specifications of Negative Pressure Control Valves

01 完美结合	电控负压控制阀与控制电位器完美的组合安装，不仅节省了安装空间，还解决了各自单独安装所带来的各种不便。	01 Perfect Combination	Between Vacuum selection valve and Potentiometer to ensure easier assembly and space saving.
02 延时功能	当工作完毕将吸唾手柄放回挂架，控制电位器即时关闭设备，3秒钟后再次启动，5秒钟后自动关闭。以便排空管内残留污物，杜绝堵塞隐患。	02 Time-delay Protection	After finished operation and put suction handle into rack, potentiometer would turn off central suction unit immediately, and restart to work for 5 seconds after 3 seconds' rest to ensure wastes inside suction tubes of vacuum system could be pumped out completely.
03 无负压启动	避免因抽吸管路负压值不够导致控制阀无法正常启动，完全解决了控制阀需要负压动能辅助启动的弊端。	03 Start without vacuum pressure	No need vacuum pressure to activate, and avoid failed activation due to lower suction power.
04 全面隔绝	内部工作腔与污物腔全面隔绝，可以有效消除产品质量隐患，并保持长期稳定工作。	04 Complete Isolation	Perfect structure to isolate wastes with main control parts to max. reduce quality risks.
05 安装灵便	根据实际需要，侧箱、地箱均可正常安装使用，无需顾虑因安装位置不当导致的质量隐患。	05 Easy Assembly	Assembly with lateral box and ground box, and no headache for place choices.
06 快速装卸	配有快速安装座，无需借助其它工具便可快速装卸。	06 Quick Dis-assembly	Equip with quick assembly and dis-assembly brackets, and no need any tools to assist.
07 适应性强	采用交/直流电源输入和抗干扰优化系统，满足现有不同中央负压抽吸设备的各种使用条件。	07 Strong Adaptability	Suitable for all kinds of central suction unit.


负压控制阀 Vacuum Suction Valve

**J38-0139**

140.6 x 101 AC/DC 24V

配置清单 Configuration List


- 电控式负压控制阀 Electronic Place Selection Valve J38-001 x 1
- 控制电位器 Potentiometer J39-001 x 1
- 接头套装 Connectors Kit J32-301 x 1
- 接头套装 Connectors Kit J32-206~J32-227 x 1
- (任选其一 Optional Parts)

**J38-001**

134.3 x 94.5 AC/DC 24V

配置清单 Configuration List

- 电控式负压控制阀 Electronic Place Selection Valve J38-001 x 1
- 接头套装 Connectors Kit J32-301 x 1
- 接头套装 Connectors Kit J32-206~J32-227 x 1
- (任选其一 Optional Parts)

**J38-010**

142.3 x 101.5

配置清单 Configuration List

- 气控式负压控制阀 Pneumatic Type Place Selection Valve J38-010 x 1
- 接头套装 Connectors Kit J32-301 x 1
- 接头套装 Connectors Kit J32-206~J32-227 x 1
- (任选其一 Optional Parts)

产品型录 PRODUCT CATALOG

单位UNIT: mm

选位阀Place Selection Valve



J06-0107

125.5 x 125.5 AC/DC 24V

配置清单Configuration List

- 电控式选位阀Electronic Place Selection Valve
- 控制电位器Potentiometer
- 接头套装Connectors Kit
- (任选其一Optional Parts)

J06-001 x 1
J07-001 x 1
J32-106~J32-127 x 1




J06-001

105 x 105 AC/DC 24V

配置清单Configuration List

- 电控式选位阀Electronic Place Selection Valve
- 接头套装Connectors Kit
- (任选其一Optional Parts)

J06-001 x 1
J32-106~J32-127 x 1



J25-001


105 x 105

配置清单Configuration List

- 气控式选位阀Pneumatic Type Place Selection Valve
- 接头套装Connectors Kit
- (任选其一Optional Parts)


J25-001 x 1
J32-106~J32-127 x 1

控制电位器 Potentiometer



J39-001

66.3 x 45.5 x 32.6 AC/DC 24V



J07-001

44.5 x 57 x 30.8 AC/DC 24V

牙科抽吸管 Dental Suction Pipes

50米/条 M/Piece(s) 100米/箱 M/Box(ES)



GZD-105-1450

Ø10.5 x 14.5mm
意大利产 Italy

GZZ-105-1450

Ø10.5 x 14.5mm
中国产 China

GZD-175-2150

Ø17.5 x 21.5mm
意大利产 Italy

GZZ-175-2150

Ø17.5 x 21.5mm
中国产 China

材质: 意大利TPO 中国PU

性能参数

1. 高耐化学性;
2. 耐温性: 长期耐温-5°C~+60°C;
3. 压力耐受: 长期耐压0.20bar, 短期耐压0.40bar。

产品特性

1. 质轻且弯曲灵活;
2. 内外壁平滑便于清洗;
3. 轴向有尼龙纱线进行加强
增加抗撕裂性和耐真空性

Materials: Italy TPO China PU

Technical Specification

1. High Chemical Resistance.;
2. Temperature Resistance: Long-Period With -5°C~+60°C.;
3. Pressure Resistance: Long-Period With 0.20bar.
Short-Period With 0.40bar.


Product Characteristics

1. Light Weight And Flexibility.
2. Easy Clean With Smooth Inner And Outer Surface.
3. Use Nylon Yarn To Reinforce Pipe On Axial Direction. Increase Tear Resistance And Vacuum Resistance.

产品型录 PRODUCT CATALOG

单位UNIT: mm

直接头Straight Connector



专注牙科设备常用配件研发制造

Focus on R & D and manufacturing of common dental equipment accessories

内接头Internal Connector



J08-006-0101



J08-008-0101



J08-010-0101



J08-012-0101



J08-013-0101



J08-014-0101



J08-016-0101



J08-018-0101



J08-019-0101



J08-020-0101



J08-021-0101



J08-022-0101



J08-025-0101



J08-027-0101




J08-110-0101

产品型录 PRODUCT CATALOG

单位UNIT: mm


管夹Pipe Clamp

100件/包 Pieces/Package



管箍Shroud Ring

50件/包 Pieces/Package




专注牙科设备常用配件研发制造

Focus on R & D and manufacturing of common dental equipment accessories


管套Pipe Shroud

100件/包 Pieces/Package





J12-122-1200

Ø2.2 x 1.2




J12-222-1200

Ø2.2 x 1.2




J12-130-1400

Ø3.0 x 1.4




J12-230-1400

Ø3.0 x 1.4





J12-140-2100

Ø4.0 x 2.1




J12-240-2100

Ø4.0 x 2.1




J12-160-4000

Ø6.0 x 4.0



J12-260-4000

Ø6.0 x 4.0




产品型录 PRODUCT CATALOG

单位UNIT: mm

接嘴Nipple


100件/包 Pieces/Package



J35-105-3300

直接嘴


Straight connection nipple



J35-905-3300

90°接嘴


connection nipple



J35-300-3200

三通接嘴

3-way connection nipple




J35-300-4832

三通接嘴


3-way connection nipple

卡扣U-Clip

O型密封圈 O Seal Ring



J05-022-0603



O型密封圈

d17.0 * 2.65mm
d18.0 * 1.5mm
d21.2 * 1.8mm
d21.2 * 2.0mm
d21.2 * 3.1mm
d28.0 * 2.0mm
d33.5 * 2.0mm

专注牙科设备常用配件研发制造

Focus on R & D and manufacturing of common dental equipment accessories

接头套装Connectors Kit



Ø6 ~ 27

J32-106

配置清单Configuration List

- 直接头Straight Connector J01-006-0101 x 1
- 内接头Internal Connector J08-006-0101 x 1
- 卡扣U-Clip J05-022-0603 x 2
- O型密封圈 O Seal Ring d21.2*1.8mm x 2



Ø6 ~ 27

J32-108

配置清单Configuration List

- 直接头Straight Connector J01-008-0101 x 1
- 内接头Internal Connector J08-008-0101 x 1
- 卡扣U-Clip J05-022-0603 x 2
- O型密封圈 O Seal Ring d21.2*1.8mm x 2



J32-110

配置清单Configuration List

- 直接头Straight Connector J01-010-0101 x 1
- 内接头Internal Connector J08-010-0101 x 1
- 卡扣U-Clip J05-022-0603 x 2
- O型密封圈 O Seal Ring d21.2*1.8mm x 2

J32-112

配置清单Configuration List

- 直接头Straight Connector J01-012-0101 x 1
- 内接头Internal Connector J08-012-0101 x 1
- 卡扣U-Clip J05-022-0603 x 2
- O型密封圈 O Seal Ring d21.2*1.8mm x 2

J32-113

配置清单Configuration List

- 直接头Straight Connector J01-013-0101 x 1
- 内接头Internal Connector J08-013-0101 x 1
- 卡扣U-Clip J05-022-0603 x 2
- O型密封圈 O Seal Ring d21.2*1.8mm x 2

J32-114

配置清单Configuration List

- 直接头Straight Connector J01-014-0101 x 1
- 内接头Internal Connector J08-014-0101 x 1
- 卡扣U-Clip J05-022-0603 x 2
- O型密封圈 O Seal Ring d21.2*1.8mm x 2

J32-116

配置清单Configuration List

- 直接头Straight Connector J01-016-0101 x 1
- 内接头Internal Connector J08-016-0101 x 1
- 卡扣U-Clip J05-022-0603 x 2
- O型密封圈 O Seal Ring d21.2*1.8mm x 2

J32-118

配置清单Configuration List

- 直接头Straight Connector J01-018-0101 x 1
- 内接头Internal Connector J08-018-0101 x 1
- 卡扣U-Clip J05-022-0603 x 2
- O型密封圈 O Seal Ring d21.2*1.8mm x 2

J32-119

配置清单Configuration List

- 直接头Straight Connector J01-019-0101 x 1
- 内接头Internal Connector J08-019-0101 x 1
- 卡扣U-Clip J05-022-0603 x 2
- O型密封圈 O Seal Ring d21.2*1.8mm x 2

J32-120

配置清单Configuration List

- 直接头Straight Connector J01-020-0101 x 1
- 内接头Internal Connector J08-020-0101 x 1
- 卡扣U-Clip J05-022-0603 x 2
- O型密封圈 O Seal Ring d21.2*1.8mm x 2

J32-121

配置清单Configuration List

- 直接头Straight Connector J01-021-0101 x 1
- 内接头Internal Connector J08-021-0101 x 1
- 卡扣U-Clip J05-022-0603 x 2
- O型密封圈 O Seal Ring d21.2*1.8mm x 2

J32-122

配置清单Configuration List

- 直接头Straight Connector J01-022-0101 x 1
- 内接头Internal Connector J08-022-0101 x 1
- 卡扣U-Clip J05-022-0603 x 2
- O型密封圈 O Seal Ring d21.2*1.8mm x 2

J32-125

配置清单Configuration List

- 直接头Straight Connector J01-025-0101 x 1
- 内接头Internal Connector J08-025-0101 x 1
- 卡扣U-Clip J05-022-0603 x 2
- O型密封圈 O Seal Ring d21.2*1.8mm x 2

J32-127

配置清单Configuration List

- 直接头Straight Connector J01-027-0101 x 1
- 内接头Internal Connector J08-027-0101 x 1
- 卡扣U-Clip J05-022-0603 x 2
- O型密封圈 O Seal Ring d21.2*1.8mm x 2

产品型录 PRODUCT CATALOG

单位UNIT: mm

接头套装Connectors Kit



Ø6 ~ 27

Ø6 ~ 27

J32-206

配置清单Configuration List

- 转接头Adapter J31-001-0101 x 1
- 直接头Straight Connector J01-006-0101 x 1
- 内接头Internal Connector J08-006-0101 x 1
- 卡扣U-Clip J05-022-0603 x 2
- O型密封圈 O Seal Ring d21.2*1.8mm x 2
- O型密封圈 O Seal Ring d21.2*3.1mm x 1

J32-208

配置清单Configuration List

- 转接头Adapter J31-001-0101 x 1
- 直接头Straight Connector J01-008-0101 x 1
- 内接头Internal Connector J08-008-0101 x 1
- 卡扣U-Clip J05-022-0603 x 2
- O型密封圈 O Seal Ring d21.2*1.8mm x 2
- O型密封圈 O Seal Ring d21.2*3.1mm x 1

J32-210

配置清单Configuration List

- 转接头Adapter J31-001-0101 x 1
- 直接头Straight Connector J01-010-0101 x 1
- 内接头Internal Connector J08-010-0101 x 1
- 卡扣U-Clip J05-022-0603 x 2
- O型密封圈 O Seal Ring d21.2*1.8mm x 2
- O型密封圈 O Seal Ring d21.2*3.1mm x 1

J32-212

配置清单Configuration List

- 转接头Adapter J31-001-0101 x 1
- 直接头Straight Connector J01-012-0101 x 1
- 内接头Internal Connector J08-012-0101 x 1
- 卡扣U-Clip J05-022-0603 x 2
- O型密封圈 O Seal Ring d21.2*1.8mm x 2
- O型密封圈 O Seal Ring d21.2*3.1mm x 1

J32-213

配置清单Configuration List

- 转接头Adapter J31-001-0101 x 1
- 直接头Straight Connector J01-013-0101 x 1
- 内接头Internal Connector J08-013-0101 x 1
- 卡扣U-Clip J05-022-0603 x 2
- O型密封圈 O Seal Ring d21.2*1.8mm x 2
- O型密封圈 O Seal Ring d21.2*3.1mm x 1

J32-214

配置清单Configuration List

- 转接头Adapter J31-001-0101 x 1
- 直接头Straight Connector J01-014-0101 x 1
- 内接头Internal Connector J08-014-0101 x 1
- 卡扣U-Clip J05-022-0603 x 2
- O型密封圈 O Seal Ring d21.2*1.8mm x 2
- O型密封圈 O Seal Ring d21.2*3.1mm x 1

J32-216

配置清单Configuration List

- 转接头Adapter J31-001-0101 x 1
- 直接头Straight Connector J01-016-0101 x 1
- 内接头Internal Connector J08-016-0101 x 1
- 卡扣U-Clip J05-022-0603 x 2
- O型密封圈 O Seal Ring d21.2*1.8mm x 2
- O型密封圈 O Seal Ring d21.2*3.1mm x 1

J32-218

配置清单Configuration List

- 转接头Adapter J31-001-0101 x 1
- 直接头Straight Connector J01-018-0101 x 1
- 内接头Internal Connector J08-018-0101 x 1
- 卡扣U-Clip J05-022-0603 x 2
- O型密封圈 O Seal Ring d21.2*1.8mm x 2
- O型密封圈 O Seal Ring d21.2*3.1mm x 1

J32-219

配置清单Configuration List

- 转接头Adapter J31-001-0101 x 1
- 直接头Straight Connector J01-019-0101 x 1
- 内接头Internal Connector J08-019-0101 x 1
- 卡扣U-Clip J05-022-0603 x 2
- O型密封圈 O Seal Ring d21.2*1.8mm x 2
- O型密封圈 O Seal Ring d21.2*3.1mm x 1

J32-220

配置清单Configuration List

- 转接头Adapter J31-001-0101 x 1
- 直接头Straight Connector J01-020-0101 x 1
- 内接头Internal Connector J08-020-0101 x 1
- 卡扣U-Clip J05-022-0603 x 2
- O型密封圈 O Seal Ring d21.2*1.8mm x 2
- O型密封圈 O Seal Ring d21.2*3.1mm x 1

J32-221

配置清单Configuration List

- 转接头Adapter J31-001-0101 x 1
- 直接头Straight Connector J01-021-0101 x 1
- 内接头Internal Connector J08-021-0101 x 1
- 卡扣U-Clip J05-022-0603 x 2
- O型密封圈 O Seal Ring d21.2*1.8mm x 2
- O型密封圈 O Seal Ring d21.2*3.1mm x 1

J32-222

配置清单Configuration List

- 转接头Adapter J31-001-0101 x 1
- 直接头Straight Connector J01-022-0101 x 1
- 内接头Internal Connector J08-022-0101 x 1
- 卡扣U-Clip J05-022-0603 x 2
- O型密封圈 O Seal Ring d21.2*1.8mm x 2
- O型密封圈 O Seal Ring d21.2*3.1mm x 1

J32-225

配置清单Configuration List

- 转接头Adapter J31-001-0101 x 1
- 直接头Straight Connector J01-025-0101 x 1
- 内接头Internal Connector J08-025-0101 x 1
- 卡扣U-Clip J05-022-0603 x 2
- O型密封圈 O Seal Ring d21.2*1.8mm x 2
- O型密封圈 O Seal Ring d21.2*3.1mm x 1

J32-227


配置清单Configuration List

- 转接头Adapter J31-001-0101 x 1
- 直接头Straight Connector J01-027-0101 x 1
- 内接头Internal Connector J08-027-0101 x 1
- 卡扣U-Clip J05-022-0603 x 2
- O型密封圈 O Seal Ring d21.2*1.8mm x 2
- O型密封圈 O Seal Ring d21.2*3.1mm x 1

专注牙科设备常用配件研发制造

Focus on R & D and manufacturing of common dental equipment accessories

接头套装Connectors Kit




J32-301

配置清单Configuration List

- 管套Pipe Shroud J12-140-2100 x 5
J12-240-2100 x 5
- 管套Pipe Shroud J12-160-4000 x 2
J12-260-4000 x 2
- 三通接嘴3-way connection nipple J35-300-3200 x 1
J35-300-4832 x 1

强力吸唾头High Volume Ejector

可调式吸唾头Adjustable Suction Connector



Ø16

J22-016


强力Strong



Ø16.5

J22-001

弱力Weak



Ø22.5


J22-002

强力Strong

产品型录 PRODUCT CATALOG


单位UNIT: mm

吸唾接头Suction Connector




J23-165

Ø6.5



J23-185

Ø8.5



J23-195

Ø9.5




J23-111

Ø11




J23-112

Ø12




J23-215

Ø15




J23-218

Ø18




J36-165

Ø6.5




J36-185

Ø8.5




J36-195

Ø9.5




J36-211

Ø11



J36-212

Ø12



J36-215

Ø15

专注牙科设备常用配件研发制造

Focus on R & D and manufacturing of common dental equipment accessories

可调式吸唾头Adjustable Suction Connector



J22-023-0165

弱力Weak




J22-023-0185

弱力Weak



J22-023-0195

弱力Weak



J22-023-0111

弱力Weak



J22-023-0112

弱力Weak



J22-023-0215

强力Strong




J22-023-0218

强力Strong

产品型录 PRODUCT CATALOG

单位UNIT: mm

可调式吸唾头Adjustable Suction Connector



J22-036-0165

弱力Weak



J22-036-0185

弱力Weak



J22-036-0195

弱力Weak



J22-036-0211

强力Strong



J22-036-0212


强力Strong




J22-036-0215

强力Strong

吸唾转换头 Suction Adapter



可调式吸唾头套装 Adjustable Suction Connector Kit




产品型录 PRODUCT CATALOG

单位UNIT: mm


吸管接头Suction Connector




J04-001-0301



J04-011-0301



J04-012-0301




J04-018-0301




J33-065

Ø6.5




J33-085

Ø8.5




J33-095

Ø9.5




J33-011

Ø11




J33-012

Ø12



J33-015

Ø15



J33-018

Ø18

直通接头 Straight Joint


45°接头 Connector



J16-001-0101



J16-002-0101





J45-025-0101

排气接头 Air Outlet Connector

90°接头 Connector



J18-006-0101



J90-025-0101

转接头 Adapter



J31-001-0101



J31-003-0101



3/4 (六分) 外管螺纹


产品型录 PRODUCT CATALOG

单位UNIT: mm


单向阀Check Valve



J10-001




J10-002



过滤网接头Filter Net Connector



J14-011



J14-021




过滤网Filter



J14-012



J14-022




过滤网Filter


专注牙科设备常用配件研发制造

Focus on R & D and manufacturing of common dental equipment accessories



三通接头T-Junction




J15-001-0101



J15-002-0101




J15-003-0101



三通接头T-Junction

封盖 Seal Cap



J15-004-0101

Ø 16* Ø 16* Ø 16mm



J15-005-0101




J17-001-0101

产品型录
PRODUCT CATALOG

单位UNIT: mm


多接头Multi-Connectors



J30-001-0101




抽吸释压阀 Vacuum Release Valve



J11-011


黑色Black



J11-012


灰色Grey

截流开关 Air Flow Adjust Valve




J09-001

Ø 6



J09-004

Ø 4



J09-006

Ø 6

专注牙科设备常用配件研发制造


Focus on R & D and manufacturing of common dental equipment accessories

Detailed specifications


详细规格

装饰管套 Decoration Casing Pipe

回气瓶 Air Return Bottle




J02-006-0201



J20-001

排污管接头 Drain Pipe Connector

*重要目的：避免医疗废气进入医疗诊室内
Main Purpose: Avoid medical disposable gas come into clinics.



J21-025-3200

Ø 25~Ø 32



Ø 25


PVC材质
PVC Material

负压抽吸专用
Only for Suction Negative Pressure



Ø 32

橡塑材质
Rubber And Plastic



J21-050-6300

Ø 50~Ø 63



Ø 50

PVC材质
PVC Material

污水排放专用
Only for Draining Negative Pressure



Ø 63

橡塑材质
Rubber And Plastic

产品型录 PRODUCT CATALOG

单位UNIT: mm

水过滤器 Water Filter



J13-011



铜烧结滤网
Copper Sintered Filter



J13-012



尼龙滤网
Nylon Filter



J13-021



铜烧结滤网
Copper Sintered Filter



J13-022



尼龙滤网
Nylon Filter

水过滤网 Water Filter Net



J13-001

尼龙滤网
Nylon Filter



J13-002

尼龙滤网
Nylon Filter



J13-003

铜烧结滤网
Copper Sintered Filter




J13-004

铜烧结滤网
Copper Sintered Filter



专注牙科设备常用配件研发制造

Focus on R & D and manufacturing of common dental equipment accessories


过滤筒Filter Cartridge



J03-001-0101





J03-011-0101





下部管接头方向通孔
Bottom connector with through-hole

过滤筒Filter Cartridge

过滤器面板Filter Panel



J03-018-0101




J03-017-0701

产品型录
PRODUCT CATALOG


单位UNIT: mm

过滤器接头Filter Connector




J03-002-0302

黄色Yellow




J03-008-0801

灰色Grey



J03-003-0302


黄色Yellow



J03-009-0302


黄色Yellow

过滤器接头Filter Connector




J03-021-0301

灰色Grey




J03-021-0302

黄色Yellow




J03-022-0301

灰色Grey




J03-022-0302

黄色Yellow



J03-019-0301

灰色Grey



J03-019-0302


黄色Yellow




产品型录 PRODUCT CATALOG

单位UNIT: mm

过滤网Filter Net




J03-012-0300




J03-020-0300


过滤网杆Disposable Filter Rod




J03-006-0302




J03-004-0302



J03-007-0302



J03-015-0302




J03-016-0302

专注牙科设备常用配件研发制造

Focus on R & D and manufacturing of common dental equipment accessories


过滤网组合套 Disposable Filter Net



J03-013-001

配置清单 Configuration List


- 过滤网 Filter Net J03-012-0300 x 1
- 过滤网杆 Disposable Filter Rod J03-006-0302 x 1



J03-014-001

配置清单 Configuration List

- 过滤网 Filter Net J03-012-0300 x 1
- 过滤网杆 Disposable Filter Rod J03-015-0302 x 1



J03-015-001


配置清单 Configuration List

- 过滤网 Filter Net J03-012-0300 x 1
- 过滤网杆 Disposable Filter Rod J03-016-0302 x 1

产品目录 PRODUCT CATALOG

单位UNIT: mm



过滤器组合套Filter Sets



J03-001

配置清单Configuration List


- 过滤筒Filter Cartridge J03-001-0101 x 1
- 过滤器接头Filter Connector J03-002-0302 x 1
- 过滤网组合套Disposable Filter Net J03-014-001 x 1
- O型密封圈 O Seal Ring d33.5*2.0mm x 1
- d18*1.5mm x 2
- d21.2*1.8mm x 1



J03-002

配置清单Configuration List


- 过滤筒Filter Cartridge J03-001-0101 x 1
- 过滤器接头Filter Connector J03-002-0302 x 1
- 过滤网组合套Disposable Filter Net J03-013-001 x 1
- O型密封圈 O Seal Ring d33.5*2.0mm x 1
- d18*1.5mm x 1
- d21.2*1.8mm x 1



专注牙科设备常用配件研发制造


Focus on R & D and manufacturing of common dental equipment accessories

过滤器组合套Filter Sets

**J03-003**


配置清单Configuration List

- 过滤筒Filter Cartridge J03-011-0101 x 1
- 过滤器接头Filter Connector J03-009-0302 x 1
- 过滤网组合套Disposable Filter Net J03-013-001 x 1
- O型密封圈 O Seal Ring d32.5*3.55mm x 1
- d21.2*1.8mm x 2

**J03-004**

配置清单Configuration List


- 过滤筒Filter Cartridge J03-011-0101 x 1
- 过滤器接头Filter Connector J03-008-0801 x 1
- 过滤网组合套Disposable Filter Net J03-014-001 x 1
- O型密封圈 O Seal Ring d21.2*1.8mm x 2



产品型录 PRODUCT CATALOG

单位UNIT: mm



过滤器组合套Filter Sets



J03-021 黄色Yellow

配置清单Configuration List



- 过滤筒Filter Cartridge J03-001-0101 x 1
- 过滤器接头Filter Connector J03-012-0302 x 1
- 过滤网组合套Disposable Filter Net J03-014-001 x 1
- O型密封圈 O Seal Ring
 - d33.5*2.0mm x 1
 - d18*1.5mm x 2
 - d21.2*1.8mm x 1




J03-031 黄色Yellow

配置清单Configuration List

- 过滤筒Filter Cartridge J03-001-0101 x 1
- 过滤器接头Filter Connector J03-022-0302 x 1
- 过滤网组合套Disposable Filter Net J03-013-001 x 1
- O型密封圈 O Seal Ring
 - d33.5*2.0mm x 1
 - d18*1.5mm x 1
 - d21.2*1.8mm x 1





过滤器组合套Filter Sets

**J03-022**

灰色Grey

配置清单Configuration List


- 过滤器筒Filter Cartridge J03-001-0101 x 1
- 过滤器接头Filter Connector J03-021-0301 x 1
- 过滤网组合套Disposable Filter Net J03-014-001 x 1
- O型密封圈 O Seal Ring d33.5*2.0mm x 1
- d18*1.5mm x 2
- d21.2*1.8mm x 1

**J03-032**

灰色Grey

配置清单Configuration List


- 过滤器筒Filter Cartridge J03-001-0101 x 1
- 过滤器接头Filter Connector J03-022-0301 x 1
- 过滤网组合套Disposable Filter Net J03-013-001 x 1
- O型密封圈 O Seal Ring d33.5*2.0mm x 1
- d18*1.5mm x 1
- d21.2*1.8mm x 1



产品类型 PRODUCT CATALOG

单位UNIT: mm


过滤器组合套Filter Sets



J03-041 黄色Yellow

配置清单Configuration List


- | | | |
|-------------------------|--------------|-----|
| • 过滤筒Filter Cartridge | J03-018-0101 | x 1 |
| • 过滤器接头Filter Connector | J03-019-0302 | x 1 |
| • 过滤网Filter Net | J03-020-0300 | x 1 |
| • O型密封圈 O Seal Ring | d28.0*2.0mm | x 1 |
| | d18*1.5mm | x 1 |
| | d21.2*1.8mm | x 1 |




J03-042 灰色Grey

配置清单Configuration List

- | | | |
|-------------------------|--------------|-----|
| • 过滤筒Filter Cartridge | J03-018-0101 | x 1 |
| • 过滤器接头Filter Connector | J03-019-0301 | x 1 |
| • 过滤网Filter Net | J03-020-0300 | x 1 |
| • O型密封圈 O Seal Ring | d28.0*2.0mm | x 1 |
| | d18*1.5mm | x 1 |
| | d21.2*1.8mm | x 1 |




过滤器组合套Filter Sets

**J03-043****黄色Yellow**


配置清单Configuration List

- 过滤器组合套Filter Sets J03-041 x 2
- 过滤器面板Filter Panel J03-017-0701 x 1
- 十字沉头螺丝Cross-shaped Screw ST2.9*6.5mm x 4

**J03-044****灰色Grey**

配置清单Configuration List


- 过滤器组合套Filter Sets J03-042 x 2
- 过滤器面板Filter Panel J03-017-0701 x 1
- 十字沉头螺丝Cross-shaped Screw ST2.9*6.5mm x 4




产品型录 PRODUCT CATALOG

单位UNIT: mm

过滤器组合套Filter Sets




过滤器组合套Filter Sets



产品型录 PRODUCT CATALOG

单位UNIT: mm

过滤器组合套Filter Sets





专注牙科设备常用配件研发制造
Focus on R & D and manufacturing
of common dental equipment accessories

Suzhou ABC Electronics Technology Co.,Ltd

Address: 7th Floor, No. 158, Kunlunshan Road, High-tech Zone, Suzhou city, Jiangsu Province, China 215 153

Website: www.abcdental.cn

Contact: Mr. Alan

Mobile: +86 152 5027 1489

Email: alan@abcdental.cn